



530



531



532



533



534



535



536



537



538



539



540



541



530.Ca



531.Ca



532.Ca



533.Ca



534.Ca



535.Ca



536.Ca



537.Ca



538.Ca



539.Ca



540.Ca



541.Ca



530	50	H	Dark green on light blue	BBa
531	80	H	Dark green on light blue	BBb
532	120	H	Dark green on light blue	BBc
533	150	H	Dark green on light blue	BBd
534	175	H	Dark green on light blue	BBe
535	225	H	Dark green on light blue	BBf
536	340	H	Dark green on light blue	BBg
537	510	H	Dark green on light blue	BBh
538	680	H	Dark green on light blue	BBi
539	750	H	Dark green on light blue	BBj
540	850	H	Dark green on light blue	BBk
541	1/2	D	Dark green on light blue	BBl
530.Ca	50	H	Inv. print d. green on light blue	
531.Ca	80	H	Inv. print d. green on light blue	
532.Ca	120	H	Inv. print d. green on light blue	
533.Ca	150	H	Inv. print d. green on light blue	
534.Ca	175	H	Inv. print d. green on light blue	
535.Ca	225	H	Inv. print d. green on light blue	
536.Ca	340	H	Inv. print d. green on light blue	
537.Ca	510	H	Inv. print d. green on light blue	
538.Ca	680	H	Inv. print d. green on light blue	
539.Ca	750	H	Inv. print d. green on light blue	
540.Ca	850	H	Inv. print d. green on light blue	
541.Ca	1/2	D	Inv. print d. green on light blue	

Set 530-541A [357 – 368]

(12 values)

Sheet of 12 of 530-541

Sheet of 12 of 530.Ca-541.Ca

Sheet of 12 of 530.Ca-541.Ca, o/p 'SPECIMEN'